



Large loop for opening or hanging.

Adhesive area on back for hanging up the package.

Instant activation - sterile salt solution included, press and it's ready to use.

Insertion guide. A firm grip without needing to touch the catheter with your fingers.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.¹



Foldable and discreet.

Smooth eyelets².

Product description

Single use, hydrophilic urinary catheter with salt solution for activation.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

Packaging

30 individually packed catheters per customer box.

Transport package: Includes either 4 or 2 customer boxes.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sterile solution pocket: Laminate foil of polyethylene terephthalate (PET), polyethylene (PE) and aluminium. Sterile water with sodium chloride (NaCl).

Item package: Polyethylene (PE).

Product range

Reference	Type	Length	Size
410 08 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 08
410 10 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 10
410 12 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 12
410 14 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 14
410 16 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 16
410 18 XX*	LoFric® Primo™ Nelaton (male)	40 cm	CH 18
411 06 XX*	LoFric® Primo™ Nelaton (Paediatric)	20 cm	CH 06
411 08 XX*	LoFric® Primo™ Nelaton (Paediatric)	20 cm	CH 08
411 10 XX*	LoFric® Primo™ Nelaton (Paediatric)	20 cm	CH 10
413 08 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 08
413 10 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 10
413 12 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 12
413 14 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 14
413 16 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 16
413 18 XX*	LoFric® Primo™ Nelaton (female)	20 cm	CH 18
414 08 XX*	LoFric® Primo™ Nelaton (female)	15 cm	CH 08
414 10 XX*	LoFric® Primo™ Nelaton (female)	15 cm	CH 10
414 12 XX*	LoFric® Primo™ Nelaton (female)	15 cm	CH 12
414 14 XX*	LoFric® Primo™ Nelaton (female)	15 cm	CH 14
415 10 XX*	LoFric® Primo™ Tiemann (male)	40 cm	CH 10
415 12 XX*	LoFric® Primo™ Tiemann (male)	40 cm	CH 12
415 14 XX*	LoFric® Primo™ Tiemann (male)	40 cm	CH 14
415 16 XX*	LoFric® Primo™ Tiemann (male)	40 cm	CH 16
415 18 XX*	LoFric® Primo™ Tiemann (male)	40 cm	CH 18

*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

LoFric® Primo™

1. Waller et al 1997; Dunger et al. 1988
2. Lundgren et al. 2000; Newman et al. 2018.

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

Wellspect. A real difference.

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

